

Balanza Analítica FA Series



LED display



Unit conversion



Multi-device printing
Standard RS232 interface



Full electromagnetic
interference protection

TECHNICAL DATA

Model	FA2204 Cal Externa	FA2204N Cal Interna
Max Capacity	200g	200g
Min Readout	0.01g	
Accuracy	0.1mg	
Repeatability	±0.0002g	
Linearity Error	±0.0003g	
Stability Time	≤ 3S	
Working Temp.	15°C ~ 35°C	
Pan Size	Standard model φ80mm	
Dimensions	300*200*295 mm	
Calibration Mode	External calibration.	/ Internal calibration
Communication Interface	RS232	

* FA-N series Internal calibration balance