

Mini Magnetic Stirrer



Strong magnetic field



Large mixing ability



Durable brushless motor



LCD display

TECHNICAL DATA

Model	MS3	MS5	MS5S
LCD Display	No	No	Yes
Control	Knob	Knob	Knob
Motor Type	Brushless motor		
Working Plate Size (mm)	125	135	135
Speed Range	350-1800rpm	350-1800rpm	200-2000pm
Max Stirring Volume	1L(H ₂ O)	2L(H ₂ O)	2L(H ₂ O)
Stir Bar Size (mm)	30	30	30
Power	5W	10W	10W
Shell Protection Level	IP21	IP42	IP42
Power adapter	DC12V 1A		